

FIXTURE REMOVER KIT




1/ FR kit

- FR kit is a tool for removing stopped implant caused by excessive torque during insertion or other several reasons.


2/ Key Points of Neo Fixture Remover kit

- When marginal bone loss occurs, implant placement and GBR can be done right after the implant removal.
- Implant can be simply removed even in broken or osseointegrated cases.
- No bone damage and no dental trauma in patients occur.
- Trephine drill is unnecessary.
- Applied to any kinds of implants.
- When implant is placed in the wrong direction, it is possible to pull implant out simply.


Fixture Remover Kit Technique Guide




Failed implant




Right turn Close




Left turn Close



Clear



Fixture Remover Kit Clinical Case




Ordering information

Catalog Number	Description
FRKIT	Fixture Remover Kit

Note

FR kit is tool for removing stopped implant causes by excessive torque during insertion or implant in lost bone easily, it's designed for implant with the same diameter after removing implant, inserting new.



Neobiotech Co.,Ltd

104, 24 / 102, 103, 104-1, 104-2, 105, 106, 205, 212, 312, 509, 510, 511, 10F,
36 Digital-ro 27-gil, Guro-gu, Seoul, Korea Tel 82-2-582-2885, Fax 82-2-582-2883
www.neobiotech.com / Made in Korea

Dongbang Acuprime

1 Forrest Units, Hennock Road East, Marsh Barton, Exeter EX2 8RU, UK
Tel +44 1392-829500 Fax +44 1392-829232

**Neo
Biotech**
Satisfaction to Dentists